TRAININGSPLAN

TRAINING PLAN

EXERCICES

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TRAININGSPLAN

EDZÉSTERV

TRÉNINKOVÝ PLÁN

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PLAN DE TRAINING

AS SEEN ON
TV
ORIGINAL

DE

Einleitung.....	03
Oszillation & Akupressur Massage.....	04
Sanfte Massage & Wadenoszillation.....	06
Tiefere Massage & Beinvibration.....	08
Wadenstärkung.....	10
Wadenmassage.....	12
Mittelfuss-Akupressur.....	14
Hintere Oberschenkelmuskulatur.....	16
IT-Bandmassage.....	18
Quadriceps.....	20
Unterer Rücken.....	22

EN

Introduction.....	03
Oscillation & Acupressure Massage.....	04
Gentle massage & calf oscillation.....	06
Deeper Massage & Leg Vibration.....	08
Calf Strengthening.....	10
Calf Massage.....	12
Metatarsal Acupressure.....	14
Rear thigh muscles.....	16
Iliotibial Band Massage.....	18
Quadriceps.....	20
Lower Back.....	22

FR

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EN INTRODUCTION

All applications should be for a maximum of 10 minutes, or as specified. You may start to feel discomfort during exercise; only use VibroLegs for as long as is comfortable to you.

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EN Oscillation & Acupressure Massage

- Sit down comfortably, lean your back against a chair or armrest, and place your feet on the foot pads.
- Use the device for a maximum of 10 minutes at a speed which is comfortable for you. This exercise is designed to relax the feet and legs, promoting circulation and muscle tension.

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EN Gentle massage & calf oscillation

- Move to the edge of the chair or sofa and sit upright while tensing your trunk. Push the appliance away far enough so that you can still place your feet completely on the pads and your calves are stretched out.
- Use the device for a maximum of 10 minutes at a speed which is comfortable for you. This exercise is less intensive for your feet and thighs, as it focuses primarily on the calves.

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EN Calf Strengthening

- Sit down and place the product as shown in Figure #1. Move the balls of your foot to the back of the pads (#2).
- Lay your toes flat on the VibroLegs and then lift your heels while tensing your calves.
- Hold this position for 5 seconds (#3), then relax for 5 seconds and lower your heels. Always keep your body taut.
- Do 3 repetitions of 10 // 15 // 20, depending on whether you are a beginner, advanced or experienced.

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EN Calf Massage

- Sit on the floor and place your calves on the pads.
- Use the VibroLegs for a maximum of 10 minutes at a speed which is comfortable for you, depending on the strength of the massage you want.

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EN Metatarsal Acupressure

- Take a comfortable sitting position and place your balls of foot between the second and third toes on the acupressure heads. Use the lowest speed and moderate pressure, so that it is comfortable.

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EN Iliotibial Band Massage

- Enter the planking position, keep your torso taut and put the VibroLegs under the side of one of your thighs, which will directly massage the iliotibial bands.
- Use the product for a maximum of 10 minutes or until you can no longer correctly hold the planking position. Use a speed that is comfortable and in proportion to the strength of the massage you want.
- Then repeat the exercise on the iliotibial band for the other leg.

NOTE: Do not place your entire body weight on the appliance, or lie down on it.

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EN Quadriceps

- Enter the planking position, keeping your torso taut, and place the VibroLegs under your thighs.
- This will directly massage the quadriceps.
- Use the product for a maximum of 10 minutes or until you can no longer correctly hold the planking position. Use a speed that is comfortable and in proportion to the strength of the massage you want. **NOTE: Do not place your entire body weight on the appliance, or lie down on it.**

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CH: Mediashop Schweiz AG | Leuholz 14 | 8855 Wangen | Switzerland
EU: MediaShop GmbH | Schneiderstraße 1, Top 1 | 2620 Neunkirchen | Austria